

# TECHNICAL TIPS

## 13

## STANDING FRAME ANGLE ADJUSTMENT

To facilitate the transfer in or out of the standing frame or to optimize the comfort of the user, the angle can be adjusted on this device. When making any changes to the angle, be sure to:

- Lock all the wheels before.
- Hold the device in place firmly, with one hand on the pivot handle while activating the cylinder activation lever with the other hand.
- Verify that no part of the user's body or that of the caregiver can be wedged or touched during the angle-change procedure.
- Check that all Velcro straps and other fasteners are in place and properly adjusted.
- Make sure the standing frame is at a proper distance away from any surrounding

To change the angle, hold the pivot handle firmly with one hand while activating the cylinder activation lever with your other hand (**fig. 13a**). To maintain positioning at a specific angle, simply release the activation lever at the desired angle. To maintain the user position at a certain angle, loosen the knob just above the cylinder activation lever and slide the round-tipped rod to the desired position to keep the angle in place (**fig. 13b**). Retighten the knob.

NOTE - If you have difficulty in reaching it or tightening it securely, use the tightening tool provided for this purpose (**fig. 13c-d**).




fig. 13a




fig. 13b




fig. 13c



fig. 13d