

S1 S2 S3

Company:

PO number:

Marked for:

BILL TO:

SHIP TO:

Address

City, province, postal code

Phone

Email



1 To ensure the Vertiflex is assembled optimally and tailored to the child's individual needs, please complete the measurement section below — **especially the 6 essential measurements highlighted**. Providing this information allows us to deliver a product that requires minimal adjustment on-site, ensuring a smoother experience for all.

Weight (pounds)

Mandatory

C Floor to axilla

E Floor to mid-hip

K Depth of axilla


L Depth of hip

H Axilla to axilla

I Hip width

Comments:

Anthropometric data



Please note that measurements C, E, H, I, K and L are essential to ensure the best experience with a Synetik stander.

2 Starter kit includes:

- Base
- Frame
- Mast
- Axilla support receiver
- Hip support receiver
- Knob adjustment tool

- Hybrid starter kit - \$ 1910 (019-00083-BH-S)
- Prone starter kit - \$ 1685 (019-00083-BP-S)
- Supine starter kit - \$ 1789 (019-00083-BS-S)

- Hybrid starter kit - \$ 1910 (019-00083-KH-S)
- Prone starter kit - \$ 1685 (019-00083-KP-S)
- Supine starter kit - \$ 1789 (019-00083-KS-S)

- Hybrid starter kit - \$ 1910 (019-00083-JH-S)
- Prone starter kit - \$ 1685 (019-00083-JP-S)
- Supine starter kit - \$ 1789 (019-00083-JS-S)

		SIZE 1			SIZE 2			SIZE 3		
		Max weight: 50 pounds			Max weight: 70 pounds			Max weight: 150 pounds		
		Standard		XL	Standard		XL	Standard		XL
		Minimum	Maximum	Maximum	Minimum	Maximum	Maximum	Minimum	Maximum	Maximum
HEIGHT	(C) Floor to axilla	17"	27"	30"	23"	33"	38"	31"	45"	53"
	(E) Floor to mid hip	10"	15"	18"	15"	20"	25"	22"	30"	38"
	(C-E) Mid hip to axilla	7"	12"	12"	8"	13"	13"	9"	15"	15"

- 3** Colour
- Purple (019-00088-PMA)
 - Turquoise (019-00088-PTU)
 - Blue (019-00088-PBL)
 - Red (019-00088-PRO)
 - Black (019-00088-PNO)

- 4** Mandatory
- Locking caster \$ 238 (019-00068)
 - Pneumatic tilt system \$ 786 (019-00069)
 - Planar pad \$ 394 (019-00038)

S1 S2 S3

5

Headrest



<p>■ 022-00005</p>	<p>Multi-angle adjustable headrest (Height x Width)</p> <p>022-00005-14 (4 ½" x 7") 022-00005-16 (5" x 8") 022-00005-17 (6" x 9")</p>	<p>\$ 638</p>
<p>■ 022-00004</p>	<p>Soft density headrest with shapeable back plate (Height x Width)</p> <p>022-00004-4 (4" x 9 ½") 022-00004-5 (5" x 10") 022-00004-6 (6" x 11")</p>	<p>\$ 350</p>
<p>■ 022-00002</p>	<p>Rigid density headrest with non-shapeable back plate (Height x Width)</p> <p>022-00002-4 (3" x 9 ½") 022-00002-5 (4" x 10") 022-00002-6 (5 ½" x 11")</p>	<p>\$ 350</p>
<p>■ 022-00001</p>	<p>Soft density headrest with non-shapeable back plate (Height x Width)</p> <p>022-00001-4 (4" x 9 ½") 022-00001-5 (5" x 10") 022-00001-6 (6" x 11")</p>	<p>\$ 350</p>
<p>■ 022-00003</p>	<p>Soft density headrest - shapeable laterals (Height x Width) 3" x 12"</p>	<p>\$ 422</p>

6

Headrest mount



<p>■ 020-00004</p>	<p>Multi-axis headrest mount 45°</p>	<p>\$ 329</p>
<p>■ 020-00001</p>	<p>Straight multi-axis headrest mount</p>	<p>\$ 329</p>
<p>■ 020-00006</p>	<p>Flip-back multi-axis headrest mount</p>	<p>\$ 467</p>
<p>■ 020-00007</p>	<p>Standard 90° headrest mount</p>	<p>\$ 467</p>
<p>■ 020-00008</p>	<p>Triple pivot multi-axis headrest mount</p>	<p>\$ 467</p>
<p>■ 020-00015</p>	<p>Mounting bracket kit - without headrest mount</p>	<p>\$ 218</p>

7


Shoulder



<p>■ 019-00006</p>	<p>Shoulder support (mandatory on supine & hybrid model)</p>	<p>\$ 202</p>
--------------------	--	---------------

8

Tray




<p>■ 019-00011</p> <p>■ 019-00010</p>	<p>Tray with half tray padding</p> <p>019-00011-20 Lexan (Clear) 019-00011-22 Lexan (Clear) 019-00011-24 Lexan (Clear)</p> <p>019-00010-22 ABS (Black) 019-00010-22 ABS (Black) 019-00010-24 ABS (Black)</p>	<p>\$ 809</p>
<p>■ 024-00010</p>	<p>Tray with half tray padding for Vertiflex mobile</p> <p>024-00010-22 Lexan (Clear) 024-00010-24 Lexan (Clear)</p>	<p>\$ 809</p>
<p>■ 019-00049</p> <p>■ 019-00057</p>	<p>Tray with activity bowl (Length 5,5"x Width 11,5"x Depth 2,25") and half tray padding</p> <p>019-00049-20 Lexan (Clear) 019-00049-22 Lexan (Clear) 019-00049-24 Lexan (Clear)</p> <p>019-00057-20 ABS (Black) 019-00057-22 ABS (Black) 019-00057-24 ABS (Black)</p>	<p>\$ 937</p>
<p>■ 019-00065</p>	<p>Elbow stop, mounted on tray (pair)</p>	<p>\$ 273</p>

S1 S2 S3

9


Upper body



019-00034	Standard thoracic support	\$ 207
019-00042	Standard thoracic support with elbow stops	\$ 717
019-00041	Swing away thoracic support	\$ 512
019-00028	Padded thoracic strap	\$ 118
019-00066	Independent thoracic lateral receiver	\$ 346
019-00044	Chest harness 019-00044-S Small 019-00044-M Medium 019-00044-L Large	\$ 315

10


Hip



019-00043	Standard hip laterals	\$ 242
019-00045	Swing away hip laterals	\$ 744
019-00035	Rigid padded hip strap (recommended for prone)	\$ 127
019-00029	Soft padded hip strap (recommended for supine and hybrid)	\$ 127
019-00008	Abduction pommel 019-00008-SS Small 019-00008-SM Medium 019-00008-SL Large	\$ 513

11


Knee



019-00033	Standard knee support with single lateral	\$ 617
019-00031	Multi-axis knee support with single lateral	\$ 739
019-00073	Standard knee support with double lateral	\$ 752
019-00071	Multi-axis knee support with double lateral	\$ 888

12

Foot



019-00084	One-piece foot positioning	\$ 504
024-00012	One-piece foot positioning for Vertiflex mobile	\$ 504
019-00085	Multi-axis foot positioning	\$ 676
140-00001	Foot holder (pair) (Width x Depth) 140-00001-1 140-00001-2 140-00001-3 140-00001-4 2 1/2" x 5 1/2" 3" x 6 3/4" 3 1/2" x 7 1/2" 3 3/4" x 8 1/2" 140-00001-5 140-00001-6 140-00001-7 140-00001-8 4" x 9 1/2" 4 1/4" x 10 1/2" 4 3/4" x 12" 5" x 13"	\$ 194
140-00001-(x)TS	Toe strap for foot holder (pair)	\$ 48

13

Other accessories



019-00002	Push handle (supine & hybrid models)	\$ 69
019-00052	Hybrid padded chinrest	\$ 350

S1 S2 S3

14

NON-REFUNDABLE
OPTION



Side wheel kit Size 2* (024-00002) Size 3* (024-00003)

1775 \$

***Refer to section 2 of this form for size selection.**

For the Vertiflex Mobile, please ensure that the N measurement (wrist to floor) is entered in section 1.